CLEAR ALIGNER

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CA® CLEAR ALIGNER AND
GINGIVAL AND PERIODONTAL
TISSUES

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All cases were carried out by Dr. Pablo Echarri. All the laboratory work was carried out by Ladent laboratory, and the Scheu Dental Technology materials were used.

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© 2017 Centro de Ortodoncia y ATM, Ladent, SL All rights reserved

This book or any part thereof may not be reproduced, stored in retrieval system or transmitted in any form or by any means electronic, mechanical, photocopying or otherwise, without prior written permission of the publishers.

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Fig 1. External root resorption grades based on Levander y cols.

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CA CLEAR ALIGNER® technique is usually prescribed to adults with periodontal alterations. The gingival and periodontal evaluation with pre-treatment panoramic x-ray is very important. Its guide in some aspects of treatment with CA CLEAR ALIGNER®.

It is necessary to have a Case History, if there some previous trauma, study the periodontum and bone support and the apical shape of the teeth or possible previous resorption.


Fig 2. Morfología radicular: Triangular (T), Romboidal (R), Pipeta (P), Dilacerada (D).

Fig 2. Root Morhology: Triangular (T), Rhomboid (R), Pipette-shaped (P), Dilacerated (D).

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Levander et al define in their study each grade of root resorption observed in radiography (Fig. 1); 0: no evidence for resorption, 1: irregular root contour, 2: apical root resorption less than 2 mm, 3: apical root resorption >2mm and <1/3 of original root length, 4: root resorption exceeding 1/3 of original root length.

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It is also important the types of root that we can observe base on its morphology (Fig. 2): triangular, rhomboid or diamond-shaped, pippete-shaped and dilacerated.

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Fig 3. Root Morhology: Triangular (T), Rhomboid (R), Pipette-shaped (P), Dilacerated (D).

It is not recommended to perform movements in more than one plan in the space in the same set-up step if patient have a reduced periodontum support. The studies demonstrate that the slight and intermittent forces cause less hyalinization zones, less root resorptions and less loss of support. So CA CLEAR ALIGNER® is ideal for those patients.

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Fig 4. Normal length of CA CLEAR ALIGNER® compared to the length indicated in patients with gingival inflammation.

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Orthodontics is the unique dentistry speciality that works with inflammation of periodontal tissues to cause an orthodontic movement. This process has to be with a controlled inflammation of the gingival tissues.

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That is why is important the hygienic control and plaque control and no try not to pressure the swollen gingival tissues with CA CLEAR ALIGNER®. If there are any inflammation the orthodontic movement has to be stopped, or alleviate the soft tissues or cut the aligner at the level of the gingival margins.

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DIASTEMA CLOSURE WITH CA
CLEAR ALIGNER®

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Diastema closure is the treatment that more gingival inflammation produces. This is why is necessary to alleviate the interdental space and soft tissues in buccal and lingual during the intrusion movement.

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Fig 6. Application of the Blue Blokker® in the to-do movement direction.

When anterior teeth have reduced periodontum support and present different axis, the patient can exerts excessive forces during insertion and removal of the CA CLEAR ALIGNER®. For this reason it has to alleviate those teeth to avoid overloads, but the alleviates have to be in the same direction as the movement.

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Fig 8. Panoramic x-ray; (A): initial, (B): final, ©: follow-up at months.

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The tooth intrusion combined with periodontal treatment and orthodontic intrusion could be a good method to improve the periodontal conditions if there are a good hygiene and force system.


Fig 9. Control de torque

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Fig 9. Torque control

The torque has to be controlled before and during the intrusion. It is not recommended intrude the tooth before a correct torque to avoid a possible root resorption.

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**INCREASING OCCLUSAL VERTICAL
DIMENSION WITH CA CLEAR
ALIGNER®**

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Fig 10. Collapsed bite because of VD loss and lingual protrusion

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The diastemas due to occlusal vertical dimension loss and collapsed bite with lingual protrusion, proinclination of upper and lower incisors and reduction of interincisal angle (Fig 10.) has to be corrected with increase of vertical dimension and retroclination of the anterior teeth.

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Fig 12. CA CLEAR ALIGNER® with extension and bubbles for lingual re-education.

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You can aggregate posterior build-ups to a CA CLEAR ALIGNER® to increase the occlusal vertical dimension (Fig 11.) and extends the upper CA CLEAR ALIGNER® to the palate making some bubbles with the CA Tip Plier 1 as the reference for the lingual re-education (Fig 12.)

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Fig 13. Summary of the treatment and case-finishing.

To increase the occlusal vertical dimension it is necessary to replace the absent teeth and restore the damaged teeth in addition to re-educate the lingual position to a more posterior zone (Fig 13.)

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PARTICULARITIES OF CA CLEAR ALIGNER® PROTOCOL WITH PERIODONTAL PATIENTS


Fig 14. A- Protocolo CA CLEAR ALIGNER® en pacientes periodontales, B- Protocolo CA CLEAR ALIGNER® standard.

Fig 14. A- Protocol of CA CLEAR ALIGNER® in periodontal patients, B- Standard protocol of CA CLEAR ALIGNER®.

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In periodontal patients you can use a variation of protocol of CA CLEAR ALIGNER® (Fig 14.)

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It is not recommended the round-treatment in periodontal patients, nor align and then make stripping. So the Progressive Stripping Technique is recommended.

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Fig 15. Scheme of supracrestal fibrotomy and the resulted periodontal fibre disposition.

Some authors recomend the supracrestal fibrotomy in major rotations that are not solved.

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DOUBLE RETENTION IN CA CLEAR
ALIGNER® TECHNIQUE

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Fig 16 . Scheme of double retention in CA CLEAR ALIGNER.

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Patients with reduced periodontum support need double retention: Fixed retention (Memotain®) and CA Retention for a longer period of time because of the tendency to relapse.

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BIBLIOGRAPHY

Melsen B, Agerbaek N, Markenstam G. Intrusion of incisors in adult patients with marginal bone loss. Am J Orthod Dentofacial Orthop 1989;96: 232-241

Melsen B, Agerback N, Eriksen J, Terp S. New attachment through periodontal treatment and orthodontic intrusion. Am J Orthod Dentofacial Orthop 1988;94:104-116

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Han JY. A comparative study of combined periodontal and orthodontic treatment with fixed appliances and Clear Aligner in patients with periodontitis. J Periodontal Implant Sci 2015;45(6)193-204

Cuoghi OA, Aiello CA, Consolaro A, Tondelli PM, Mendonça MR de. Resorption of roots of different dimension induced by different types of forces. Braz Oral Res. 2014;28(1):1-7.

Levander E, Malmgren O. Evaluation of the risk of root resorption during orthodontic treatment: a study of upper incisors. Eur J Orthod 1988;10:30-38.