CLEAR ALIGNER

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All cases were carried out by Dr. Pablo Echarri. All the laboratory work was carried out by Ladent laboratory, and the Scheu Dental Technology materials were used.

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© 2018 Centro de Ortodoncia y ATM, Ladent, SL All rights reserved

This book or any part thereof may not be reproduced, stored in retrieval system or transmitted in any form or by any means electronic, mechanical, photocopying or otherwise, without prior written permission of the publishers.

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Consider:

Facial pattern.
"Oral" display.
Occlusal plane rotation.

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Treatment protocol

Forces.
Torque control: before and during the intrusion.
Movement sequence.

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Fig. 1. Types of gingival insertion; (A): Normal, (B): Normal with deep pockets, and (C): Gingival insertion at an enamel level.

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A normal gingival insertion (Fig. 1A) doesn't require treatment, normal gingival insertion with deep pocket (Fig. 1B) is treated with gingivectomy and the insertion at the enamel level (Fig. 1C) requires crown lengthening with apical repositioning flap.

In front of vertical problem like anterior deep bite, we treat dolicofacial patients with anterior teeth intrusion. Braquifacial patients are treated with molar extrusion (Fig. 2).

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Fig. 2. (A): Molar extrusion in braquifacial patients, and (B): incisal intrusion in dolicofacial patients.

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We must bear in mind that, in both situations, treatment in the upper arch provokes an counterclockwise rotation of the occlusal plane and the treatment in lower arch provokes a clockwise rotation of the occlusal plane (Fig. 3).

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Fig. 3. Oclusal plane rotation according to the arch treated.

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We also work with light forces to avoid the blood supply cessation (Fig. 4) and correct the torque before the intrusion to ensure the root movement through the cancellous bone (Fig. 5).

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Fig. 5. Radicular problems due to no-control of the torque before the intrusion.

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It will be necessary the use of Power grips 8 and 9 (Figs. 6 and 7) during the intrusion to avoid the torque increasing of the upper incisors. Using the power grips we apply the intrusion force in the same direction than the tooth's axis (Fig. 8).

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Figs. 6 and 7. On the left, CA Power grip 8 - BUAS (Build-up Anterior Small) for lateral incisors. On the right, CA Power grip 9 - BUAL (Build-up Anterior Large) for central incisors.

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Fig. 8. Comparison between an intrusion movement with and without power grips.

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The placement of CA Power grip requires the CA Holder (Fig. 10) and corresponding mold (Figs 11 and 12) based on the scheme of Fig. 14.

Fig. 10. CA Holder

Fig. 10. Ca Holder.

Fig. 11. CA Power grip 8 - BUAS.


Fig. 11. CA Power grip 8 - BUAS.

Fig. 12. CA Power grip 9 - BUAL.

Fig. 12. CA Power grip 9 - BUAL.

Fig. 13. CA Power grip 8 colocado en el CA Holder.

Fig. 13. CA Power grip 8 into the CA Holder.

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Fig. 14. Scheme of measures of CA Power grip 8 and 9.

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
PLACEMENT PROTOCOL FOR CA
POWER GRIP 8 and 9

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1- Placement of CA Power Grip at 2 mm from the uncial ridge.

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

We produce the "splint effect" using CA Power grip 8 and 9, stimulating the closing musculature contraction.

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USING CA CLEAR ALIGNER WITH CA
POWER GRIP 8 and 9

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Clear Aligner before the intrusion. It doesn't fit properly in the molars.

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Clear Aligner for intrusion. After the intrusion it fits correctly.

Fig. 16. Fotografías iniciales.

Fig. 16. Inicial photos.


Fig. 17. Registros digitales iniciales.

Fig. 17. Initial digital registration.

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Fig. 18. Progress 1 - Alignment with proinclination and stripping.

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Fig. 19. Progress 2 - Alignment with proinclination and stripping.

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