CLEAR ALIGNER

25

CA[®] CLEAR ALIGNER & PERIODONCIA

CA[®] CLEAR ALIGNER & PERIODONTICS

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All cases were carried out by Dr. Pablo Echarri. All the laboratory work was carried out by Ladent laboratory, and the Scheu Dental Technology materials were used.

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© 2016 Centro de Ortodoncia y ATM, Ladent, SL
All rights reserved

This book or any part thereof may not be reproduced, stored in retrieval system or transmitted in any form or by any means electronic, mechanical, photocopying or otherwise, without prior written permission of the publishers.

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CA[®] CLEAR ALIGNER & PERIODONTICS

Before starting any orthodontic treatment, and therefore, CA[®] CLEAR ALIGNER treatment too, it is very important to carry out the evaluation of periodontal and muco-gingival health of the patient.

It should be evaluated:

- A. If a patient presents gingival inflammation.
- B. If he presents bleeding.
- C. If he presents bacterial plaque.
- D. If he presents tartar.
- E. If he presents gingival sulcus depth more than 2 mm deep (pathological pocket).
- F. If he presents gingival recession.
- G. The thickness of soft tissues.
- H. If he presents vertical and horizontal bone resorption.

Inflammation and bleeding should be treated before starting the orthodontic treatment. Bacterial plaque and tartar should also be removed before starting the treatment. If a gingival recession is present, expansion treatment should not be carried out. The soft tissue thickness is very important, because thick gingiva biotype reacts with inflammation and overgrowth, and a thin biotype reacts with recession. In this way it is possible to predict the behavior of soft tissues.

If any changes of patient's gingival health are observed, it is important to carry out the evaluation of his periodontal health again.

All the studies indicate that the accumulation of bacterial plaque is less in patients who use removable orthodontic appliances (including aligners), than in patients who use fixed orthodontic appliances, because they have better access when brushing the teeth and less food retention.

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There is a direct relation between bacterial plaque and gingival inflammation. Because of this, it is very important to recommend correct hygiene of CA® CLEAR ALIGNER to the patient with CETRON® and to show him the correct brushing technique (modified Bass technique).

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Figure 1. Brushing with modified Bass technique.

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Figure 2. CA® CLEAR ALIGNER which usually covers 2 mm of soft tissues in order to improve the esthetics and to apply a force closer to the rotation center of the tooth, achieving in this way more effective tooth movement.

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Figure 3. Clean CA® CLEAR ALIGNER with CETRON® following the instructions.

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When 1x1 protocol is followed, one step of the treatment is carried out on the basis of one impression, but in 1x2 and 1x3 protocols, two or three steps respectively are carried out on the basis of the same impression. The first step of any of these protocols is carried out with CA® CLEAR ALIGNER splints covering 2 mm of gingiva, but in second and third step, CA® CLEAR ALIGNER splints has to be trimmed up to gingiva because it is not possible to predict with exactness the behavior of soft tissues.

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
Figure 4. CA® CLEAR ALIGNER covering 2 mm of soft tissues in 1x1 protocol, or in the first step of 1x2 and 1x3 protocol.

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Figure 5. CA® CLEAR ALIGNER trimmed up to gingival margins, used in the second step of 1x2 protocol and the second and the third step of 1x3 protocol.

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If a patient presents gingival inflammation it is better to not cover the soft tissues with CA® CLEAR ALIGNER, and trim it up to the necks or block-out it with BLUE BLOKKER® so it doesn't press the papillae.

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Figure 6. A patient with gingival inflammation.

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
Figure 7. CA® CLEAR ALIGNER trimmed up to the necks due to gingival inflammation.

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Figure 8. CA® CLEAR ALIGNER blocked-out in the necks due to gingival inflammation.

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Space closure treatments can provoke the inflammation of papillae in closed spaces. To prevent this inflammation, the correct tooth brushing is very important and it may be also necessary to carry out the prophylaxis before taking the impressions.

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Figure 9. Space closure and inflammation of papillae.

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In anterior deep bite treatment with incisors intrusion it is very important to block-out labial and lingual or palatal soft tissues with BLUE BLOKKER® before molding the CA® CLEAR ALIGNER.

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
Figure 10. Deep bite case in which the incisors and canine intrusion is indicated.

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Figure 11. Set-up model for intrusion of incisors and canines.

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
Figure 12. Labial and palatal block-out of soft tissues in incisors and canines with BLUE BLOKKER®.

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
Figure 13. Molding of CA® CLEAR ALIGNER over the model.

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Figure 14. Insertion of CA® CLEAR ALIGNER in the mouth. Observe that, at the beginning, CA® CLEAR ALIGNER is separated from occlusal surfaces of posterior teeth and anterior soft tissues. After a couple of days, the intrusion of the teeth will take place, and the appliance will complete its function. It will not carry out excessive pressure over the soft tissues because they were previously blocked-out.

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In stripping treatments it is very important to take into account the studies of Tarnow which establish that the papilla is 4.5 mm long, measured from the bone crest to the vertex of papilla.

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According to Tarnow, if interdental contact point is 5 mm away from the bone crest, the patient presents 100% of papillae. If the contact point is 6 mm away from the bone crest, the papillae will be present in only 56% of cases. If the contact point is 7 or more millimeters away from the bone crest, the papilla will be present in 25% of cases.

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Figure 16. The relation between the position of the contact point and the presence of interdental papillae.

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In this way, in patients who present black gingival triangles, it is important to establish the contact point 5 mm away from the bone crest in order to assure the presence of papillae. This can be achieved with stripping and reapproximation of the teeth in crowding and Bolton discrepancy excess cases, or with interproximal esthetic reconstruction in the cases in which the patient doesn't present Bolton discrepancy.

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Figure 17.
 A. Presence of black gingival triangles.
 B. Stripping and reapproximation.
 C. Esthetic reconstructions.

CLINICAL PAPER

CLEAR ALIGNER