CLEAR ALIGNER

26

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WHEN IS THE RIGHT MOMENT TO USE
THE CA® POWER GRIP?

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All cases were carried out by Dr. Pablo Echarri. All the laboratory work was carried out by Ladent laboratory, and the Scheu Dental Technology materials were used.

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© 2016 Centro de Ortodoncia y ATM, Ladent, SL
All rights reserved

This book or any part thereof may not be reproduced, stored in retrieval system or transmitted in any form or by any means electronic, mechanical, photocopying or otherwise, without prior written permission of the publishers.

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CA® WHEN IS THE RIGHT MOMENT TO USE THE CA® POWER GRIP?

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INTRODUCTION

This Clinical Paper explains the right timing for the use of each of CA® POWER GRIP forms.

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CA® POWER GRIP 1

The Ca® POWER GRIP 1 - MFM (Multifunctional Forced Movement) is used to correct rotations and to extrude teeth.

These movements are carried out after the space has been obtained to facilitate the tooth position correction, and therefore, the CA® POWER GRIP is bonded after the space gaining.

Rotated teeth usually don't allow the necessary access for the CA® POWER GRIP bonding before the space gaining.

The space is achieved by usual methods of CA® CLEAR ALIGNER technique: expansion or stripping, and when this stage is completed, this is the best moment for bonding the CA® POWER GRIP 1 and correction of the rotated tooth. To correct rotations, it is necessary to try to bond the MFM at the mesio-distal center of the tooth. The impression is taken with already bonded MFM to carry out CA® CLEAR ALIGNER Forced Rotation.

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Figure 1.

- A: Case with rotated incisor in which there's no access to bond el CA® POWER GRIP 1.
- B: Expansion treatment with CA® CLEAR ALIGNER to obtain the space for correction of the rotation.
- C: After the expansion, there is access to bond CA® POWER GRIP 1.
- D: CA® POWER GRIP 1 bonded.
- E: CA® CLEAR ALIGNER Forced Rotation.
- F: Correction is finished.
- G: Finished case.

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If the access space is still not sufficient to use MFM, CA[®] TIP 2 Plier can be used to correct rotation. The necessary space for the movement is blocked out during the preparation of the set-up model, and the CA[®] TIP 2 Plier is used to carry out the "activation points" which will facilitate the rotation movement.

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Figure 2.

- A: Case with rotated incisor in which there's no access to bond el CA[®] POWER GRIP 1.
- B: Expansion treatment with CA[®] CLEAR ALIGNER to obtain the space for correction of the rotation.
- C: Preparation of the model blocking out the necessary zones to carry out the correction.
- D: CA[®] CLEAR ALIGNER.
- E: Formation of "activation points" with CA[®] TIP 2 Plier.
- F: Correction with CA[®] CLEAR ALIGNER activado.
- D: Finished case.

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In extrusion cases, the correction space should be gained, too, as a previous step, and by means of expansion and stripping. Then, MFM is bonded and the access is always possible. The MFM for extrusion should be bonded in the mesio-distal center of the tooth and in a close position to gingival margin. The impression is taken with bonded MFM and CA® CLEAR ALIGNER Forced Extrusion is carried out.

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Figure 3.

- A: Case with intrusion of upper left central incisor and without space to carry out extrusion.
- B: Expansion with CA® CLEAR ALIGNER.
- C: Bonding of CA® POWER GRIP 1.
- D: CA® POWER GRIP 1 bonded.
- E: CA® CLEAR ALIGNER Forced Extrusion with elastic.
- F: Extrusion is carried out.
- G: Finished case.

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When the correction is completed, the CA® POWER GRIP 1 can be removed.

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CA[®] POWER GRIP 2 & 3

The CA[®] POWER GRIP 2 - MDI (Mesio-Distal Inclination) and the CA[®] POWER GRIP 3 - DMI (Disto-Mesial Inclination) are used to correct mesio-distal inclinations of the teeth.

The CA[®] POWER GRIP 2 facilitates the correction of the teeth with a root movement from mesial to distal in upper right and lower left quadrant. In upper left and lower right quadrants, it facilitates the root movement from distal to mesial.

The CA[®] POWER GRIP 3 facilitates the correction of the teeth with a root movement from distal to mesial in upper right and lower left quadrant. In upper left and lower right quadrants, it facilitates the root movement from mesial to distal.

In both cases, first the space must be obtained, necessary for the correction by means of expansion and/or stripping. Then the CA[®] POWER GRIP 2 or 3 is bonded in the mesio-distal center of the tooth and the impression is taken to carry out the following step of CA[®] CLEAR ALIGNER.

When the correction is completed, the CA[®] POWER GRIP 2 or 3 can be removed.

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Figure 4.

A: Case with distal inclination of upper left central incisor and without space for correction.
 B: Expansion with CA® CLEAR ALIGNER.
 C: Bonding of the CA® POWER GRIP 3.
 D: CA® POWER GRIP bonded.
 E: Correction with CA® CLEAR ALIGNER.
 F: Correction is finished.
 G: Finished case.

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CA[®] POWER GRIP 4 & 5

The CA[®] POWER GRIP 4 - FALA (Force Application Labial) and 5 - FALI (Force Application Lingual) are bonded close to the gingival margin to correct the bucco-lingual inclination or torque of the teeth. If they are bonded on labial surface of the teeth they increase the torque, and if they are bonded on lingual surface of the teeth, they reduce the torque.

If the teeth are crowded or rotated, they don't offer the necessary access for bonding of the CA[®] POWER GRIP, but anyway, the torque should not be corrected in this stage of the treatment. The first step is to gain space for correction of teeth position, and then the alignment and correction of rotations of the teeth. This is the moment for bonding the CA[®] POWER GRIP 4 or 5 on the corresponding surface of the tooth, then the impression is taken, and the CA[®] CLEAR ALIGNER is fabricated. The access for CA[®] POWER GRIP bonding is always possible. Nevertheless, some teeth with reduced size of microdontic teeth might have insufficient surface for CA[®] POWER GRIP 5 bonding, which is smaller. In such cases, CA[®] TIP Plier 4 is used to carry out "activation points" in CA[®] CLEAR ALIGNER, properly prepared for this activation (i.e., during the set-up preparation, the space for tooth correction was planned in the opposite side of the aligner activation side).

When the correction is completed, the CA[®] POWER GRIP 4 or 5 can be removed.

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Figure 5.
 A: Case with crowding and buccal position of the lower right central incisor.
 B: Expansion carried out with CA® CLEAR ALIGNER.
 C: Alignment finished with negative torque of lower right central incisor.
 D: Bonding of CA® POWER GRIP 5.
 E: CA® POWER GRIP bonded.
 F: Correction with CA® CLEAR ALIGNER.
 G: Finished case.

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Figure 6.

A: Case with crowding and buccal position of the lower right central incisor.
 B: Expansion carried out with CA® CLEAR ALIGNER.
 C: Alignment finished with negative torque of lower right central incisor.
 D: Fabrication of "activation points" with the CA® TIP 4 Plier.
 E: Correction with CA® CLEAR ALIGNER.
 F: Finished case.

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CA® POWER GRIP 6 & 7

The CA® POWER GRIP 6 - BUB (Build-Up Bicuspid) and 7 - BUM (Build-Up Molar) are used to carry out posterior build-ups with an aim to increase vertical dimension and to facilitate the correction of crossed teeth. They can be used both with CA® CLEAR ALIGNER and with fixed orthodontic appliances.

They are bonded when the position of crossed teeth is going to be corrected, and before taking the impression to carry out the CA® CLEAR ALIGNER. The CA® POWER GRIP 6 and 7 are bonded to bicuspids and molars and the occlusion should be adjusted with articulating paper in order to obtain the four contact points which would lift the occlusion in a stable way.

When the correction of cross bite is completed, the CA® POWER GRIP 6 or 7 can be removed.

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Figure 7.

A: Case with crowding and anterior cross bite.
 B: Bonding with the CA® POWER GRIP 6 and 7 in lower molars and bicuspids.
 C: CA® POWER GRIP 6 and 7 bonded.
 D: Correction with CA® CLEAR ALIGNER.
 E: Correction is finished.
 F: Finished case.

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CA[®] POWER GRIP 8 & 9

The CA[®] POWER GRIP 8 - BUAS (Build-Up Anterior Small) and 9 - BUAL (Build-Up Anterior Large) are used in cases with anterior deep bite with normal overjet. They are bonded to the palatal surface of the upper central and lateral incisors to facilitate both intrusion of incisors and extrusion of molars.

When the incisors present rotations or crowding, it may be impossible to access for bonding of the CA[®] POWER GRIP, but anyway, this is not the right moment for intrusion of the teeth. First the space must be obtained by means of expansion and/or stripping, the crowding and rotations must be corrected, and the torque must be corrected to make sure the intrusion of incisors would be carried out within the spongy bone without the interference of the cortical bone.

When the teeth are aligned, the access for CA[®] POWER GRIP bonding is sufficient. When the bonding is carried out on the palatal surface of the upper incisors and at the distance of 1.5 - 2 mm from their incisal edge, the occlusion is adjusted with the articulating paper so the lower incisors make contact with the four upper CA[®] POWER GRIP (the rest of the teeth remain out of occlusion), and the impression is taken to carry out the following step of the CA[®] CLEAR ALIGNER treatment.

These anterior build-ups should not be removed after the completion of correction. When the treatment is finished, the palatine surface form of the upper incisors is completed with composite, creating in this way the necessary cingula which will act as permanent retention. If the build-ups are removed, the lower incisors would remain out of occlusion which would provoke the movement of extrusion resulting in the relapse of deep bite.

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Figure 8.

A: Case with crowding and anterior deep bite in which it is not possible to bond CA® POWER GRIP 8 and 9.
 B: Expansion carried out with CA® CLEAR ALIGNER.
 C: Correction of crowding and rotations with CA® CLEAR ALIGNER.
 D: Bonding of CA® POWER GRIP 8 and 9.
 E: Correction of anterior deep bite.
 F: Correction is finished, keeping the anterior build-ups as permanent retention.

CLINICAL PAPER

CLEAR ALIGNER