CLEAR ALIGNER

29

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TREATMENT WITH DISTALIZATION
USING CA[®] CLEAR ALIGNER
Y CA[®] DISTALIZER KIT
(Part 1)

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All cases were carried out by Dr. Pablo Echarri. All the laboratory work was carried out by Ladent laboratory, and the Scheu Dental Technology materials were used.

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© 2016 Centro de Ortodoncia y ATM, Ladent, SL
All rights reserved

This book or any part thereof may not be reproduced, stored in retrieval system or transmitted in any form or by any means electronic, mechanical, photocopying or otherwise, without prior written permission of the publishers.

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
TREATMENT WITH DISTALIZATION
USING CA® CLEAR ALIGNER
Y CA® DISTALIZER KIT
(Part 1)

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Figs. 1-5. Case 1402. It presents mild molar Class III malocclusion with crowding in both jaws.

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Figs. 6-7. X-rays of lower molars. Observe the presence of the right lower third molar and the absence of the left lower third molar.

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Figs. 8-11. The treatment in maxilla is carried out with CA® CLEAR ALIGNER for transverse expansion with the VECTOR® 40 screw. When the space is obtained, the treatment is continued with sequential splints for alignment.

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Figs. 12-15. In mandible, CA® CLEAR ALIGNER with CA® Distalizer Kit acc. Dr. Echarri (see Ca® Clinical Paper 27 and 28) is used. The activation of the screw is carried out with the frequency of 1 turn each 2-3 days.

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
Figs. 16-19. Molar Class I is achieved after 4-week treatment.

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Figs. 20-21. Observe that the right side is more distalized than the left side. This difference is due to the presence of the bud of the third molar on the right side. The distalization is 2 mm on the right side, and 2.5 mm on the left side.

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Figs. 22. CA® Distalizer Kit with VECTOR® 350 Screw inserted in CA® CLEAR ALIGNER

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Figs. 23-26. When the molar Class I is achieved, the distalization is finished. CA® POWER GRIP forms should be bonded in canines and bicuspid, then an impression should be taken to be able to continue the treatment. CA® POWER GRIP 3 is bonded in lower right bicuspid and canine, and Ca® POWER GRIP 2 is bonded in lower left bicuspid and canine. Models with CA® POWER GRIP forms in canines and bicuspid.

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Figs. 27-29. CA CLEAR ALIGNER during the distalization of bicuspids.

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Fig. 30. Canine and bicuspid distalization is carried out simultaneously if it is not necessary to correct rotation, to align or to level these teeth. If it is necessary, they should be distalized one by one, and the malpositioned ones should be corrected

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Fig. 31. Continuation of distalization of canines and bicuspid.

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Fig. 32. When the canine Class I is achieved, the incisors alignment is carried out with progressive set-ups. In continuation, the final detailing is carried out.

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Figs. 33-36. Models after the distalization of molars and with Ca® POWER GRIP 2 and 3 in canines and bicuspid.

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Fig. 37. Model on model holder.

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Figs. 38-39. Using Ca® SMART 4.0, it is possible to take occlusal, frontal, right lateral and left lateral records. To fix the position of the chess plate, fixing blocks are used. The photos show occlusal calibration and left calibration.

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
Fig. 40. Initial photograph of right side. The molars were distalized 2 mm and 2 bicuspids will be distalized simultaneously in 2 steps (1 mm per step), and then the canine will be distalized in 2 steps. The distalization of the left side is carried out in the same way.

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
Figs. 41-43. Model on the model holder with horizontal cut.

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Figs. 44-45. First set-up with distalization of right bicuspids. They were distalized 1 mm and the angulation was overcorrected 5°.

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Fig. 46. Record of the first set-up of the right side, distalizing the bicuspids.

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
Fig. 47. Measurement of distalization of 1 mm.

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Fig. 48. Measurement of angulation overcorrection of 5°.

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Figs. 49-50. First set-up with distalization of left bicuspids. They were distalized 1 mm and the angulation was overcorrected 5°.

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Fig. 51. Record of the first set-up of left side distalizing the bicuspids.

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
Fig. 52. Measurement of distalization of 1 mm.

Fig. 53. Measurement of angulation overcorrection of 5°.

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WITH THE DISTALIZATION OF 1 mm OF LEFT BICUSPIDS AND RIGHT BICUSPIDS, THE ALIGNERS ARE CARRIED OUT.

**CONTINÚA EN
CA® CLINICAL PAPER 30**

**CONTINUES IN
CA® CLINICAL PAPER 30**

CLINICAL PAPER

CLEAR ALIGNER