CLINICAL PAPER

  
CLEAR ALIGNER

  
CLEAR ALIGNER

**36**

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LAST CONTRIBUTIONS TO  
CA® CLEAR ALIGNER

  
INTEGRATED SOLUTIONS IN ORTHODONTICS

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*All cases were carried out by Dr. Pablo Echarri. All the laboratory work was carried out by Ladent laboratory, and the Scheu Dental Technology materials were used.*

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© 2017 Centro de Ortodoncia y ATM, Ladent, SL All rights reserved

This book or any part thereof may not be reproduced, stored in retrieval system or transmitted in any form or by any means electronic, mechanical, photocopying or otherwise, without prior written permission of the publishers.

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The CA VECTOR® 45 has been design for expansion in the upper arch with CA CLEAR ALIGNER® hard. With the two metallic guides offer more stability to the aligner during the expansion of anterior area.

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Fig 2. CA CLEAR ALIGNER® aesthetic

It was designed hard, soft and medium foils in A2 colour, specially for patients with missing teeth or bleaching need. Due to aesthetics, it is recommended to trim at the level of gingival margins.

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### CA® STRIP OPENER

Fig 3. CA® STRIP OPENER

The CA® STRIP OPENER is a field to carry out stripping specially design for the cases when it is not possible to create a separation between the teeth before performing the stripping. There is available in 40 µm (red ring) and 15 µm (orange ring).

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The CA® CLEAR ALIGNER E-FUNCTIONAL is indicated in skeletal class II cases, with mandibular retrognatia, meso or braquifacial patients. Includes in the CA® CLEAR ALIGNER hard the ADVANCING WINGS, plastic buttons for class II elastics and it is possible to include CA® VECTOR SCREWS for expansion and in the set-up model carry out the molar extrusion.


Fig 5. Esquema del tornillo Tractive CA VECTOR® 350

Fig 5. Scheme of Tractive CA VECTOR® 350 screw

Fig 6. Tornillo Tractive VECTOR® 350.

Fig 6. Tornillo Tractive VECTOR® 350.

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**ACTIVACIÓN**

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**INDICATION**

Mesialization of upper and lower molars.

**POSITION**

One-way screw in lingual zone and tube, sectional arch in buccal zone.

**OBJETIVE**

Mesialization of upper and lower molars.

**ACTIVATION**

1/4 turn of the screw each 2 to 3 days.