CLEAR ALIGNER

16

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*All cases were carried out by Dr. Pablo Echarri. All the laboratory work was carried out by Ladent laboratory, and the Scheu Dental Technology materials were used.*

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© 2014 Centro de Ortodoncia y ATM, Ladent, SL  
All rights reserved

This book or any part thereof may not be reproduced, stored in retrieval system or transmitted in any form or by any means electronic, mechanical, photocopying or otherwise, without prior written permission of the publishers.

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## CA Clear Aligner: Rotations Correction

CA Clear Aligner Technique allows the correction of rotations using different systems.

1. CA Clear Aligner and progressive set-ups
2. CA Clear Aligner and CA Power Grip
3. CA Clear Aligner and Activation Points (CA Plier Tip 2)

### ROTATIONS CORRECTION WITH CA CLEAR ALIGNER

CA CLEAR ALIGNER & PROGRESSIVE SET-UPS

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Figure 1. Rotations Correction with CA Clear Aligner.

**CA Clear Aligner y set-ups  
progresivos**

**CA Clear Aligner and  
progressive set-ups**

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Fig. 2. Case with rotated lower right central incisor.

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Fig. 3. The space necessary for rotation correction by means of expansion and stripping has to be achieved.

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Fig. 4. Rotation correction of 0.5 mm.

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Fig. 5. Rotation correction of 1 mm.

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Fig. 6. The A part shows a scheme in which the distal surface of the central incisor is correct, but the mesial surface of the tooth should be corrected by means of its contraction. In each step of set-up, 1 mm should be corrected, until the complete correction is achieved.

The B part shows a scheme in which the distal surface of the central incisor is correct, but the mesial surface of the tooth should be corrected by means of its expansion. In each step of set-up, 1 mm should be corrected, until the complete correction is achieved.

In the C part, the rotation should be corrected by means of contraction of the distal surface of the tooth, and by means of expansion of its mesial surface. A 0.5 mm is corrected at each side.

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Fig. 7. CA Power Grip 1 (MFM Multifunctional Force Movement) is used for forced rotation.

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Fig. 8. Case with rotated lower right central incisor.

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Fig. 9. The space necessary for the rotation correction by means of expansion or stripping has to be achieved.

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Fig. 10. CA Power Grip 1 is bonded on the lower right central incisor.

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Fig. 11. A set-up is carried out with a 0.5 mm of rotation correction.

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Fig. 12. A set-up with 1 mm of rotation correction is carried out, until a complete rotation correction is achieved.

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Fig. 13. The scheme of rotation correction by means of forced rotation.

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Fig. 14. When the CA Power Grip 1 is bonded, the impressions are taken, and all set-ups necessary for complete rotation correction are carried out (see CA Newsletter 7). A set-up with 0.5 mm correction is carried out, and CA Clear Aligner SOFT and MEDIUM are molded over it (1S - 1M). Another set-up with 1 mm correction is carried out over the same model, and CA Clear Aligner SOFT and MEDIUM are molded over it (2S - 2M). The fabrication of set-ups is continued until the complete correction of rotation is achieved (3S - 3M).

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Fig. 15. Using the same initial model, the space necessary for complete correction is foreseen, and a CA Clear Aligner HARD Forced Rotation (HFR) is moulded, trimmed at the level of the CA Power Grip, and with a CA Lingual Button. The patient wil use it with an elastic from the button of the aligner to the CA Power Grip. The elastic should be of 1/8" and 2.5 oz, or 4.5 oz, and the patient should use a new elastic each time he uses the aligner.

ROTATION

		
1 <sup>st</sup> week	1 S	HFR
2 <sup>nd</sup> week	1 M	HFR
3 <sup>rd</sup> week	2 S	HFR
4 <sup>th</sup> week	2 M	HFR
5 <sup>th</sup> week	3 S	HFR
6 <sup>th</sup> week	3 M	HFR

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Fig. 16. The first week, the patient will use CA Clear Aligner 1S when he is outside, and CA Clear Aligner HFR at home and at night. The second week, he will use 1M when he is outside, and the same HFR at home and at night. He will continue in the same way with 2S, 2M, 2S, 3M aligners.

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Fig. 17. A case with rotated upper left central incisor. The space for rotation correction has already been achieved by means of stripping.

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Fig. 18. Progressive set-up for rotation correction.

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Fig. 19. Superimposition of initial model and set-up.

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Fig. 20. When the CA Clear Aligner is inserted into the mouth, it will press on one side and there will be some space on the opposite side.

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Fig. 21. When the correction is achieved, the tooth will fill in the foreseen space.

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Fig. 22. If you want to carry out an activation of an aligner, you should prepare an extra space on the set-up model to facilitate the correction.

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Fig. 23. When the aligner is inserted in the mouth, the space for correction is increased.

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Fig. 24. When CA Clear Aligner HARD has carried out its effect, there's still some extra space left.

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Fig. 26. After correction, carried out by the Activation Point.

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Fig. 27. Ativation with the CA Plier Tip 2.

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Fig. 28. Protocol used with CA Plier Tip 2. Request SOFT and MEDIUM aligners from the lab, as well as HARD aligner with some extra space for activation. The first week, a patient will use the SOFT aligner, the second week he will use the MEDIUM one, and the third week he will use the HARD. Then he has to go to the clinic to take impressions for the next step. The clinician will activate the HARD aligner using the CA Plier Tip 2, and the patient will use this activated aligner during the fourth week.

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Fig. 29. If you don-t have the plier, you can request SOFT, MEDIUM, HARD and HAP (Hard with Activation Point) aligners from the lab.

1 <sup>st</sup> week	SOFT	
2 <sup>nd</sup> week	MEDIUM	
3 <sup>rd</sup> week	HARD	
4 <sup>th</sup> week	HAP	

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Fig. 30. The first week, a patient will use the SOFT aligner, the second week, the MEDIUM, and the third week he will use the HARD, then he will go to the clinic to take impressions. During the fourth week, he will use the HAP.

**NEWSLETTER**

Dr. Pablo Echarri